

# Introducing The New KitCheck Tags



Provides a faster scanning experience.



Unlock the better functionality of next generation RFID tags.



75% smaller than legacy tags



No RFID Printer required

## Why are we sunsetting the old tags?

- The new tags are safer because they do not have printed medication information that could be out-of-sync with the drug label itself
- The new tags have been engineered to improve the speed & performance
- Most hospitals are no longer using printer tags

## New Tags vs. Legacy Tags

### New Tags

Size: 0.7x1.0 (Flag)

Workflow:

- (1) Tag meds of the same NDC/Lot/Exp
- (2) Place group of tagged meds into station
- (3) Hit the Scan button
- (4) Input med info into App

### Legacy Tags Workflow

Size: 1.1x2.5 (Flag)

Workflow:

- (1) Input NDC in app
- (2) Enter in the LOT/EXP into app
- (3) Input in the quantity
- (4) Print tags & Tag Items